





Drawn by:	KA	Date	26-11-2025
Revised	- /KA	Date	-
Revised	- /KA	Date	-
Copyright reserved by law		Unit of measurement: mm	

Overview parts Reflex XXL

Lopital Nederland bv

Draw. no: EXPV5900-2025-01



P:\SW-Tek\Reflex XXL 5900\

Drawn by:	KA	Date	26-11-2025
Revised	- /KA	Date	-
Revised	- /KA	Date	-

Copyright reserved
by law

Unit of measure-
ment: mm

Overview parts Reflex XXL

Lopital Nederland bv

Draw. no: EXPV5900-2025-01

Reflex XXL, Article no: 5100 5900

Part no.	Art.no.	QTY.	Name	Drawing no.
1	70590001	1	Chassis	5900-11
2	70590002	1	Seat frame	5900-12
3	70590003	1	Main lever	5900-13
4	70590004	2	Control lever	5900-14
5	70590005	2	Seat support	5900-15.1
6	70590012	1	Battery Cover	5900-20.1
7	70590035	2	Calf strap holder	5900-24
8	70590013	12	Threaded bush Inox Ø12xM8x63 mm	5900-21
9	70590020	4	Threaded bush Inox Ø12xM8x23 mm	5900-29
10	70590029	2	Threaded bush Inox Ø12xM8x29 mm	5900-34
11	70590037	1	Bracket controlbox	5900-26
12	70590032	1	Bracket handcontrolplug	5900-36
13	70590015	1	Guide for po set XXL	5900-25
14	70590034	2	Cable duct	5900-35
15		1	Head support	Head support XL XXL
16		1	Armrest set width adjustment XL/XXL	Armlegger breedteverstelling XL XXL
17		1	Footrest set width adjustment XL/XXL	Voetsteunset XXL
18	70590009	2	Seat part XXL	5900-30
19	70590010	1	Backrest XXL	5900-31
20	70560066	4	Protective cover castor	5600-48
21	70519931	2	Clamp set footrest blue	5500-24
22	70590014	2	Cover Seat Support	5900-22
23	70590019	1	Calf Strap XXL	5900-28
24	70590022	1	Direction castor double roll Ø125 braked	Blickle Zwenkwiel XXL
25	70590021	3	Castor double roll Ø125 braked	Blickle Zwenkwiel XXL
26	70595022	1	Side Cover set	5900-33
27	70200041	2	Pump/actuator cap	3D_TA23-Kap
28	70590036	1	Control box	LOP-2022 behuizing + beugel Assy V2024-04-28
29	70590023	2	Actuator LA23 Tilt	LA23 - Reflex XL XXL
30	70590024	2	Actuator LA40 High/Low	Linak LA40 HL
31	70550018	1	Chassispart charger	Lader Magn
32	70560004	1	Battery Set	Accu 12V 7Ah
33	70590030	1	Handcontrol	3D_TMH3 TBV Reflex XL XXL
34	70590038	1	Emergency stop	YW1B-V4E02R EMERGENCY STOP (IDEC)
35		1	Sticker E-Stop	5900-94
36	70590025	1	Control Panel	Bedieningspaneel
37	70550012	1	Charger	Oplader Magn
38	70590033	1	Cable set XL XXL	Kabelset XL XXL
39	70560040	20	Bearing bush brass Ø16xØ12x4 mm	5600-56
40	70560041	20	Bearing bush PA Ø24,15xØ12,1x23 mm	5600-57.5
41	70602006	8	Snaploc	Snaploc 4029 007 0601
42	70820019	4	Cross Knob PA 40 mm M6 grey	Kruisknop M6
43	70590049	3	Wing screw M8	Kipp K0274.0108X20

Drawn by: KA	Date 26-11-2025
Revised - /KA	Date -
Revised - /KA	Date -

Copyright reserved
by law

Unit of measure:
mm


Overview parts Reflex XXL

Lopital Nederland bv

Draw. no: EXPV5900-2025-01

Reflex XXL, Article no: 5100 5900

Part no.	Art.no.	QTY.	Name	Drawing no.
44		2	Tube cap Ø30	<i>Insteekdop rond 32 grijs</i>
45		2	Flexicap round PVC 8x15 grey	<i>Dompeldop Rond PVC 8 X 15 Grijs</i>
46		12	Screw M5x6 ISO7380 A2	<i>M3x6 ISO7380</i>
47		2	Screw M4x12 DIN7985 A2	<i>M4x12 DIN7985</i>
48		4	Screw M5x6 ISO7380 A2	<i>M5x6 ISO7380</i>
49		4	Screw M5x8 ISO7380 A2	<i>M5x8 ISO7380</i>
50		4	Screw M5x10 ISO7380 A2	<i>M5x10 ISO7380</i>
51		8	Screw M6x10 DIN965 A2	<i>M6x10 DIN965</i>
52	70590044	36	Screw M8x10 ISO7380 + Dry-loc+Lurinox	<i>M8x10 ISO7380-2 Dry-loc+Lurinox</i>
53		4	Screw M8x35 DIN912 A2	<i>M8x35 DIN912</i>
54		4	Screw M10x80 DIN912 A2	<i>M10x80 DIN912</i>
55		2	Locknut M5 DIN985	<i>M5 DIN985</i>
56		4	Washer M5 Ø15 DIN9012 A2	<i>Carrosseriering M5 DIN9021 A2</i>
57		28	Washer M8 Ø25	<i>Carrosseriering M8x20</i>
58		4	Nord-Lock M10	<i>Nord-Lock M10</i>



Drawn by:	KA	Date	26-11-2025
Revised	- /KA	Date	-
Revised	- /KA	Date	-

Copyright reserved
by law

Unit of measure-
ment: mm



Overview parts Reflex XXL

Lopital Nederland bv

Draw. no: EXPV5900-2025-01

Head support XL XXL

Part no.	Art.no.	QTY.	Name	Drawing no.
1	70590045	1	Headrest Frame	5900-41
2	70560056	1	Adjuster depth	Versteller hoofdsteen Reflex
3	70560064	1	Head/necksupport fastener	Montage hoofdsteen
4	70590031	1	PUR head/neck support	PUR hoofdsteen XXL - Froli
5	70590049	1	Wing screw M8	Kipp K0274.0108X20
6		3	Screw M5x10 ISO7380 A2	M5x10 ISO7380
7		1	Tube cap 15x15	Insert 15x15



Drawn by:	KA	Date	27-11-2025
Revised	- /KA	Date	-
Revised	- /KA	Date	-

Copyright reserved
by law

Unit of measure-
ment: mm


Overview parts Head support

Lopital Nederland bv

Draw no.: EV327-2025-01

Armrest set width adjustment XL/XXL

Part no.	Art.no.	QTY.	Name	Drawing no.
1	70590011	2	Adjustment bar armrest	5900-19.2
2	70590006	1	Frame armrest XL XXL set	5900-16.2LR
3	70590046	2	Round bar armrest	5900-37.1
4	70590047	2	Unlock button armrest	5900-38
5	70820013	2	Cap	8200-23.1
6	70570503	2	Armrest PU	Froli armsteun - 5457-50 32 mm
7		2	Washer M25 DIN433 PA	M25 DIN 433 PA
8	70590048	2	Spring 0,8x11x80 Inox	Veer 0,8x11x80
9		2	Tube cap Ø28	Insteedkop rond 28 grijs
10		2	Circlip 25 DIN471	25 DIN 471
11		2	Screw M6x10 DIN914 A2	M6x10 DIN 914
12		2	Screw M5x10 DIN912 A2	M5x10 DIN912
13		2	Screw M6x16 ISO7380 A2	M6x16 ISO7380
14		4	Screw M6x35 DIN7991 A2	M6x35 DIN7991
15		2	Cylindric pin Ø8x36 DIN7	8x36 DIN 7 m6
16		2	Cylindric pin Ø8x30 DIN7	8x30 DIN 7



Drawn by: KA	Date 27-11-2025
Revised - /KA	Date -
Revised - /KA	Date -

Copyright reserved
by law

Unit of measure:
mm




Overview parts Armrest set width adjustment XL/XXL

Lopital Nederland bv

Draw
no.: EXPVArmleggerBVXLXXL-2025-01

Footrest set width adjustment XL/XXL

Part no.	Art.no.	QTY.	Name	Drawing no.
1	70590028	2	Frame width adjustment footrest	5900-18
2	70590027	2	Frame footrest XXL	5900-17
3	70590040	2	Frame footplate XXL	5900-40
4	70530032	1	Footplate set	5300-37LR.2
5	70530029	2	Pull knob M12x1.5	Treknop M12x1.5 - Kort
6		2	Tube cap 45x25	Insteekdop 25x25
7		2	Screw M5x10 DIN912 A2	M5x10 DIN912



Drawn by:	KA	Date	27-11-2025
Revised	- /KA	Date	-
Revised	- /KA	Date	-

Copyright reserved
by law

Unit of measure-
ment: mm

Overview parts
Footrest set width adjustment XL/XXL
Lopital Nederland bv

Draw.
no.: EXPVVoetsteunetXXL-2025-01